

sempermed® 
derma plus

DERMA PLUS PRODUCED SAFETY

HIGH TEAR RESISTANCE – through double-immersion production process

PERFECT GRIP – through micro-rough surface

ANATOMICALLY PERFECT SHAPE – with broader back of the hand
and naturally curved finger


- EST. 1920 -
sempermed®
YOUR GLOBAL GLOVE EXPERT™

PRODUCED SAFETY

Manufactured while using the specially developed double-dipping process the sempermed® derma plus excels with good donning properties and an AQL of 0.65.

PRODUCT HIGHLIGHTS

- High tear resistance and reliable impermeability through double-immersion production process.
- Good grip through micro-rough surface texture, even in a moist environment.
- Anatomically perfect shape with broader and more comfortable back of the hand and naturally curved fingers.



TECHNICAL DATA

TYPE	sterile surgical glove for single use, lightly powdered			
MATERIAL	natural rubber latex (NR)			
COLOUR	natural white			
GLOVE SHAPE	fully anatomical with rolled rim			
SIZE/OVERALL LENGTH as per EN 455-2	5½, 6 and 6½ / median 270mm 7, 7½ and 8 / median 280mm 8½ and 9 / median 285mm			
WALL THICKNESS measured in single layer in the palm area	0.19-0.24 mm			
BARRIER PERFORMANCE Tightness exceeding EN 455-1 (AQL 0.65 instead of AQL 1.5)	AQL 0.65			
FORCE AT BREAK exceeding EN 455-2 (≥ 12N instead of 9N)	median 12N			
DURABILITY in original packaging if stored as per DIN 7716, ISO 2230	5 years			
STERILISATION	radiation STERILE R			
PACKAGING	left and right glove with turned up cuff in inner pouch, ozone-tight, sealed in medical peel pack. in dispenser carton with sterilisation indicator: 50 pairs in transport carton with sterilisation indicator: 300 pairs			
LABELLING	EN 1041, EN 455, EN 420, EN 374, ISO 15223			
PURPOSE as set out in MD Directive 93/42/EEC as set out in PPE Regulation (EU) 2016/425	medical device class IIa single-use protective glove – PPE category III (protection against chemical substances for limited time)*			
ARTICLE NUMBERS	size 5½	823152529	size 7½	823152729
	size 6	823152609	size 8	823152809
	size 6½	823152629	size 8½	823152829
	size 7	823152709	size 9	823152909



CONTACT US!

Semperit Investments Asia Pte. Ltd.
8 Jurong Town Hall Road, #29 – 03 to 06 The JTC Summit,
Singapore 609434, Singapore
Tel.: +65 6274 4861 · Fax.: +65 6274 6977

Sempermed Europe GmbH
Am Belvedere 10,
A-1100 Vienna, Austria
Tel.: +43 1 79 777-0 · Fax: +43 1 79 777-630

sempermed@semperitgroup.com · www.sempermed.com

* Please note the Sempermed Chemical Resistance list is available on <https://www.sempermed.com/knowledge-area/chemical-resistance/chemical-resistance-sempermed-brands/>
IMPORTANT NOTE: The latest product information is available at www.sempermed.com under the «Products» tab. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. **Semperit is not liable for incorrect use of the gloves nor will Semperit be liable for improper storage and handling. In case of damage or unusual signs of wear and tear, gloves shall not be used.** In case of any doubt, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors and amendments. **CAUTION: Natural latex can cause allergic reactions, including anaphylactic shock. Semperit bears no liability for allergic reactions.** ©Copyright 2022 Sempermed