

sempermed® 
derma PF

DERMA PF FOR SENSITIVE HANDS

SKIN-FRIENDLY – no powder-related allergies
EASY TO DON – patented synthetic inner coating
WORK CONFIDENTLY – micro-rough surface for great grip



100 
+ YEARS
- EST. 1920 -
sempermed®
YOUR GLOBAL GLOVE EXPERT™

FOR SENSITIVE HANDS

The Sempermed® Derma PF offers good reasons to switch to powder free gloves.

■ SKIN-FRIENDLY

This powder free glove provides improved safety for staff, patients and visitors as there are no powder-related dermatological, allergic and toxic health risks.

■ EASY TO DON

The glove is easy to put on thanks to the synthetic lining.

■ SURE GRIP

The Sempermed® Derma PF offers a confident grip and thus combines safety and tactile properties.

■ AN ECONOMICALLY SOUND SOLUTION

Allergies and other health risks not only pose risks to end users, they can also put a significant strain on health care systems. This powder free glove removes the risk of powder-related health problems and thereby reduces potential costs in the long term.



TECHNICAL DATA

TYPE	sterile surgical glove for single use, powder free, with synthetic lining			
MATERIAL	natural rubber latex (NR)			
COLOUR	natural white			
GLOVE SHAPE	fully anatomical with rolled rim			
SIZE/OVERALL LENGTH as per EN 455-2	5½, 6 and 6½ / median 260mm 7, 7½ and 8 / median 270mm 8½ and 9 / median 280mm			
WALL THICKNESS measured in single layer in the palm area	0.16-0.21 mm			
BARRIER PERFORMANCE Tightness exceeding EN 455-1 (AQL 1.0 instead of AQL 1.5)	AQL 0.65			
FORCE AT BREAK exceeding EN 455-2 (≥ 12N instead of 9N)	median 12N			
DURABILITY in original packaging if stored as per DIN 7716, ISO 2230	3 years			
STERILISATION	radiation STERILE R			
PACKAGING	left and right glove with turned up cuff in inner pouch, ozone-tight, sealed in flat medical fibre-free peel pack. in dispenser carton with sterilisation indicator: 50 pairs in transport carton with sterilisation indicator: 300 pairs			
LABELLING	EN 1041, EN 455, EN 420, EN 374, ISO 15223			
PURPOSE as set out in MD Directive 93/42/EEC as set out in PPE Regulation (EU) 2016/425	medical device class IIa single-use protective glove – PPE category III (protection against chemical substances for limited time)*			
ARTICLE NUMBERS	size 5½	822751522	size 7½	822751722
	size 6	822751602	size 8	822751802
	size 6½	822751622	size 8½	822751822
	size 7	822751702	size 9	822751902



CONTACT US!

Semperit Investments Asia Pte. Ltd.
8 Jurong Town Hall Road, #29-03/06 The JTC Summit,
Singapore 609434, Singapore
Tel.: +65 6274 4861 · Fax.: +65 6274 6977

Sempermed Europe GmbH
Am Belvedere 10,
A-1100 Vienna, Austria
Tel.: +43 1 79 777-0 · Fax: +43 1 79 777-630

sempermed@semperitgroup.com · www.sempermed.com

* Please note the Sempermed Chemical Resistance list is available on www.sempermed.com/chemical-resistance-brands

IMPORTANT NOTE: The latest product information is available at www.sempermed.com under the «Products» tab. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. **Semperit is not liable for incorrect use of the gloves nor will Semperit be liable for improper storage and handling. In case of damage or unusual signs of wear and tear, gloves shall not be used.** In case of any doubt, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors and amendments. **CAUTION: Natural latex can cause allergic reactions, including anaphylactic shock. Semperit bears no liability for allergic reactions.** ©Copyright 2022 Sempermed